MSDS

Material Safety Data Sheet

Date prepared: August, 29th, 2012 Date revised: N/A (first edition)

Section 1 - General Information

Chemical Name: Didecyl phthalate

Synonyms: NA **CAS Reg. Number:** 84-77-5 **Catalog Number:** 1ST1132 Manufacturer's Name: A Chemtek, Inc.

100 Barber Avenue

Worcester, MA 01606, USA

Telephone Number: 508-471-4121 Fax Number: 508-845-9201

For R & D use only

Section 2 - Hazardous Ingredients/Identity Information

Classification of the substance or mixture

Not a hazardous substance or mixture

GHS Label elements, including precautionary statements

Not a hazardous substance or mixture.

Section 3 - Composition/Information on Ingredients

Formula: $C_{28}H_{46}O_4$ Molecular Weight: 446.66 g/mol CAS-No. 84-77-5

Section 4 - First Aid Measures

GENERAL ADVICE: Consult a physician. Show this safety data sheet to the doctor in attendance.

Move out of dangerous area.

EYE CONTACT: Rinse thoroughly with plenty of water for at least 15 minutes and consult a

physician.

SKIN CONTACT: Wash off with soap and plenty of water. Consult a physician.

If breathed in, move person into fresh air. If not breathing, give artificial INHALATION:

respiration. Consult a physician.

INGESTION: Never give anything by mouth to an unconscious person. Rinse mouth with

water. Consult a physician.

Section 5 - Fire-Fighting Measures

Flammable properties

Flash point 200 °C (392 °F) - closed cup

Ignition temperature not available

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

Special protective equipment for fire-fighters

Wear self-contained breathing apparatus for firefighting if necessary.

Section 6 - Accidental Release Measures

Personal precautions Use personal protective equipment. Avoid breathing vapours, mist

or gas. Ensure adequate ventilation. Evacuate personnel to safe areas.

Environmental precautions Prevent further leakage or spillage if safe to do so. Do not let

product enter drains.

Methods for cleaning up Sweep up and shovel. Keep in suitable, closed containers for

disposal.

Section 7 - Handling and Storage

Handling Avoid inhalation of vapour or mist.

Storage Keep container tightly closed in a dry and well-ventilated place.

MSDS

Material Safety Data Sheet

Date prepared: August, 29th, 2012 Date revised: N/A (first edition)

Section 8 - Exposure Controls/Personal Protection

Components with workplace control parameters

Personal protective equipment

Respiratory protection Where risk assessment shows air-purifying respirators are appropriate use a

full-face respirator with multi-purpose combination (US) or type ABEK (EN 14387) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection Handle with gloves. Gloves must be inspected prior to use. Use proper

glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and

dry hands.

Eye protection Safety glasses with side-shields conforming to EN166 Use equipment for

eye protection tested and approved under appropriate government standards

such as NIOSH (US) or EN 166(EU).

Skin protection Impervious clothing, The type of protective equipment must be selected

according to the concentration and amount of the dangerous substance at the

specific workplace.

Hygiene measures Handle in accordance with good industrial hygiene and safety practice.

Wash hands before breaks and at the end of workday.

Section 9 - Physical and Chemical Properties

Appearance: Form Liquid

Color light yellow

Safety data

pH: not available Melting point: not available

Boiling point: 265 - 268 °C (509 - 514 °F) Flash point: 200 °C (392 °F) - closed cup

Ignition temperature: not available Lower explosion limit: not available Upper explosion limit: not available Water solubility: not available

Section 10 - Stability and Reactivity

Storage stability Stable under recommended storage conditions.

Materials to avoid Strong oxidizing agents.

Hazardous decomposition products Hazardous decomposition products formed under fire conditions.

Carbon oxides.

Section 11 - Toxicological Information

Acute toxicity:

Oral LD50

LD50 Oral - rat - > 61,824 mg/kg

Inhalation LC50 not available Dermal LD50

LD50 Dermal - rabbit - 16,229 mg/kg Other information on acute toxicity

no data available

Irritation and corrosion: not available

Sensitization: not available **Respiratory or skin sensitization:** not available

MSDS

Material Safety Data Sheet

Date prepared: August, 29th, 2012 Date revised: N/A (first edition)

Chronic exposure No component of this product present at levels greater than or

equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC, ACGIH, NTP, or OSHA

Section 12 - Ecological Information

Elimination information (persistence and degradability)

not available

Ecotoxicity effects: not available

Further information on ecology: not available

Section 13 - Disposal Considerations

Product Offer surplus and non-recyclable solutions to a licensed disposal company.

Contact a licensed professional waste disposal service to dispose of this

material.

Contaminated packaging: Dispose of as unused product.

Section 14- Transport Information

DOT (US):

Not dangerous goods

IMDG:

Not dangerous goods

IATA:

Not dangerous goods

Section 15 - Regulatory Information

SARA 302 Components

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

No SARA Hazards

Massachusetts Right To Know Components

No components are subject to the Massachusetts Right to Know Act.

Pennsylvania Right To Know Components

Didecyl phthalate CAS-No. 84-77-5 Revision Date

New Jersey Right To Know Components

Didecyl phthalate CAS-No. 84-77-5 Revision Date

California Prop. 65 Components

This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

Section 16 - Other Information

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. A Chemtek shall not be held liable for any damage resulting from handling or from contact with the above product.